# Detaseal® hydroflow lite regular/fast

VPS precision impression silicone, low consistency – light bodied

Fig. 1

## 1. Tray preparation

Select a tray type suited for desired application (standard-size or custom, in metal or hard plastics). To avoid teeth contact with the tray, place stops from wax or light-curing resin into the tray at non-relevant areas of the impression. For optimal adhesion of impression material, we recommend applying a thin layer of **trayloc® A** to tray (Fig. 1). Let dry for approx. 3 minutes.

Fig. 2

## 2. Mixing and dispensing

Insert cartridge into Automix2-System mixing gun (Fig. 2). Remove and dispose of closure cap. For exact flow control extrude slightly material until uniform flow from both orifices is achieved (Fig. 3). Attach mixing cannula by aligning notched side to cartridge and turning clockwise until fastened (Fig. 4). Dispenser is now ready to mix silicone in any amount needed. After each use leave the used mixing cannula in place to seal the cartridge.


Fig. 3

## 3. Application

Apply **Detaseal® hydroflow lite** according to impression technique (e.g. double-mix technique, Fig. 5). For applications in the groove area use intra-oral-tip or disposable syringe (Fig. 6) Then seat the impression until **Detaseal® hydroflow lite** has set.


Fig. 4

## 4. Disinfecting

After removing tray from the mouth, rinse impression under running, lukewarm water. If desired, impression can then be disinfected by immersing into 2% glutaraldehyde for 15 minutes or other usual solutions.


Fig. 5

## 5. Pouring the model

We recommend waiting at least 30 minutes before pouring the model. The impression will retain its dimensional stability for up to 14 days. All class III and IV dental plasters and standard modelling resins can be used.

Fig. 6

## 6. Cleaning the tray

After setting, material can be removed with a blunt instrument. Remaining film of **trayloc® A** is dissolved by immersion in acetone or isopropanol (be sure to use only in well-ventilated rooms). Tray can then be cleaned and disinfected in the usual manner.

## Important working hints


- Substances such as retraction liquids can adversely affect the setting reaction and should therefore be thoroughly rinsed off before taking the impression.
- Do not combine with condensation curing silicones.
- Contact with latex examination gloves or other latex material may impair the setting reaction of **Detaseal® hydroflow lite**.
- Be sure to remove all traces of impression material from patient's mouth.
- Cured impression materials are chemically inert – spots on clothing should be avoided.
- We recommend use of **Detaseal® hydroflow putty, -heavy** or **-Jumbo putty** as tray material.

Further information:


Silicone based materials are proven a million times. On condition of a proper application, undesired effects are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to make an allergy test before the application of the material.

For dental use only.

### Detaseal® hydroflow lite regular

### Detaseal® hydroflow lite fast

## Indications for use:

- Impressions with sandwich and double mix technique
- Corretion impressions
- Reline impressions

## Technical data:

- DIN EN ISO 4823 - Type 3
- Mixed volume:** 50 ml cartridge
  - Mixing ratio:** 1:1
  - Product colour:** base: pink catalyst: white
  - Mixing time:** omitted (Automix2-System)
  - Working time:** **regular:** approx. 2 min. \* **fast:** approx. 1 min. \*
  - Setting time:** **regular:** approx. 4 min. \* **fast:** approx. 2 min. 30 sec. \*
  - Time in the mouth:** **regular:** approx. 2 min. **fast:** approx. 1 min. 30 sec
  - Strain in compression:** approx. 6,5 %
  - Recovery from deformation:** > 99,8 %
  - Linear dimensional change:** < 0,2 %
  - Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity.
  - Storage:**

\* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity. Increased temperatures accelerate, decreased temperatures retard these times.

## Ordering information:

**Detaseal® hydroflow lite**  
Standard packing  
2 cartridges of 50 ml  
12 mixing cannulas, yellow  
**regular set** 02730  
**fast set** 02733

**Detaseal® hydroflow Xlite**  
Standard packing  
2 cartridges of 50 ml  
12 mixing cannulas, yellow  
**regular set** 02741  
**fast set** 02744

Mixing cannulas 02706  
yellow, 48 pcs.

**Detaseal® hydroflow putty**  
Standard packing 02727  
2 jars of 250 ml  
2 measuring scoops

**Jumbo packing** 03283  
maxi cartridge of 300 ml base + 60 ml catalyst  
6 dynamic mixers, yellow  
1 rotating ring, yellow

**Detaseal® hydroflow heavy**  
Jumbo packing 02729  
maxi cartridge of 300 ml base + 60 ml catalyst  
6 dynamic mixers, yellow  
1 rotating ring, yellow

**trayloc® A** 03098  
17 ml bottle with brush

**DETAX**

GmbH & Co. KG  
Postfach 100225 · 76256 Ettlingen / Germany  
Telefon 07243/510-0 · Fax 07243/510-100  
www.detax.de · mail@detax.de

Made in Germany

IM05

02/2013

# Detaseal® hydroflow lite regular/fast

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