


Hartmetallbohrer • Tungsten Carbide Burs • Fraises en Carbure de Tungstène

H1 (001)

Rund
Round
Rond



Größe/Size/taille ø 1/10 mm	005	006	008	010	012	014	016	018	021	023	027
Länge-Length in mm - longueur											


Shank   

US. No.		1/4	1/2	1	2	3	4	5	6	7	8	10
FG 500 314 001 001...	H1- ... -FG	005	006	008	010		014		018		023	
RA 500 204 001 001...	H1- ... -RA		006	008	010	012	014	016	018	021	023	027

FG  opt. 160.000 rpm RA  opt. 5.000 - 40.000 rpm

H1S (001)

Rund, schnittfreudig
Round, fast cutting
Rond, puissant



Größe/Size/taille ø 1/10 mm	006	007	008	009	010	012	014	016	018	021	023	025	027	029
Länge-Length in mm - longueur														


Shank   

		US. No.	1/2S		1S		2S	3S	4S	5S	6S	7S	8S	9S	10S
FG		500 314 001 003...	H1S- ... -FG	006	007	008	009	010	012	014	016	018	021	023	
FGXL		500 316 001 003...	H1S- ... -FGXL					010	012	014	016	018	021	023	
RA		500 204 001 003...	H1S- ... -RA					010	012	014	016	018	021	023	025 027 029
RAL		500 205 001 003...	H1S- ... -RAL						012	014	016	018	021	023	

 opt. 160.000 rpm FG/FGXL  opt. 8.000-10.000 rpm RA/RAL

H1SX (001)

Rund, Kreuzverzahnung
Round, Cross Cut, fast cutting
Rond, Denture croisée



Größe/Size/taille ø 1/10 mm	010	012	014	016	018	021	023
Länge-Length in mm - longueur							


Shank   

US. No.		2SX	3SX	4SX	5SX	6SX	7SX	8SX
RA  500 204 001 190...	H1SX- ... -RA	010	012	014	016	018	021	023
RAL  500 205 001 190...	H1SX- ... -RAL		012	014	016	018	021	023

 opt. 8.000-10.000 rpm

H1SF (004)

Rund
Round
Rond



Größe/Size/taille ø 1/10 mm	014	016	018
Länge-Length in mm - longueur			


  

RA  500 204 004 009...	H1SF-...-RA	014	016	018
--	--------------------	------------	------------	------------

 opt. 1.500 - 2.000 rpm

H1SXN (001)

Rund, Kreuzverzahnung,
schlanker Hals
Round, Cross Cut, Slim Neck
Rond, denture croisée
avec un col qui est long et fin



Größe/Size/taille ø 1/10 mm	010	012	014	016	018	021	023
Länge-Length in mm - longueur							

Shank   

US. No.		2SXN	3SXN	4SXN	5SXN	6SXN	7SXN	8SXN
RA  500 204 001 190...	H1SXN- ... -RA	010	012	014	016	018	021	023
RAL  500 205 001 190...	H1SXN- ... -RAL			014		018		023

 opt. 1.000 - 1.500 rpm

