


Diamantscheiben • Diamond Discs • Disques diamantés

Superflex
806 104 ...

355



Größe/ Size/ Taille	ø 1/10 mm	190	220
Belegung/ Coating/ Granulation	B mm	3,00	3,00
Kopf/ Head Length/ Longueur mm	L mm	0,12	0,12


Shank   

HP **806.104. ...** **355.504.190** **355.504.220**

ø 190 ⤵ opt. 30.000 rpm • ø 220 ⤵ opt. 25.000 rpm

Superflex
806 104 ...

355



Größe/ Size/ Taille	ø 1/10 mm	160	190	220
Belegung/ Coating/ Granulation	B mm	3,00	3,00	3,00
Kopf/ Head Length/ Longueur mm	L mm	0,15	0,15	0,15



Shank   

HP **806.104. ...** **355.514.160** **355.514.190** **355.514.220**

ø 160/ø 190 ⤵ opt. 30.000 rpm • ø 220 ⤵ opt. 25.000 rpm

Superflex
806 104 ...

355



806.104.355.524.220



806.104.355.514.160



Größe/ Size/ Taille	ø 1/10 mm	190	220
Belegung/ Coating/ Granulation	B mm	3,00	3,00
Kopf/ Head Length/ Longueur mm	L mm	0,25	0,25

Shank   

HP **806.104. ...** **355.524.190** **355.524.220**

ø 190 ⤵ opt. 30.000 rpm • ø 220 ⤵ opt. 25.000 rpm

	Gelber Ring/Yellow Ring/Bague jaune	Superfein/Superfine/Superfin	30 µm	ISO No. 504
	Roter Ring/Red ring/Bague rouge	Fein/Fine/Fin	50 µm	ISO No. 514
	Blauer Ring/Blue ring/Bague bleue	Standard/Standard/Standard	100 µm	ISO No. 524