






# Diamantscheiben • Diamond Discs • Disques diamantés

Flex

806 104 ...

321



| Größe/ Size/ Taille            | 190  | 220  | 190  | 220  |
|--------------------------------|------|------|------|------|
| Belegung/ Coating/ Granulation | 3,00 | 3,00 | 3,00 | 3,00 |
| Kopf/ Head Length/ Longueur mm | 0,20 | 0,20 | 0,30 | 0,30 |

Shank   



HP **806.104. ...** ■ **321.514.190** ■ **321.514.220** ■ **321.524.190** ■ **321.524.220**

190  opt. 30.000 rpm • 220  opt. 25.000 rpm

Flex

806 104 ...

335



| Größe/ Size/ Taille            | ø 1/10 mm | 220         |
|--------------------------------|-----------|-------------|
| Belegung/ Coating/ Granulation | B mm      | voll - full |
| Kopf/ Head Length/ Longueur mm | L mm      | 0,30        |

Shank   




HP **806.104. ...** ■ **335.524.220**

 opt. 25.000 rpm

Flex

806 104 ...

345



| Größe/ Size/ Taille            | ø 1/10 mm | 190         | 220         |
|--------------------------------|-----------|-------------|-------------|
| Belegung/ Coating/ Granulation | B mm      | voll - full | voll - full |
| Kopf/ Head Length/ Longueur mm | L mm      | 0,20        | 0,20        |

Shank   

HP **806.104. ...** ■ **345.514.190** ■ **345.514.220**

190  opt. 30.000 rpm • 220  opt. 25.000 rpm


|                                                                           |                              |        |             |
|---------------------------------------------------------------------------|------------------------------|--------|-------------|
| <span style="color: yellow;">■</span> Gelber Ring/Yellow Ring/Bague jaune | Superfein/Superfine/Superfin | 30 µm  | ISO No. 504 |
| <span style="color: red;">■</span> Roter Ring/Red ring/Bague rouge        | Fein/Fine/Fin                | 50 µm  | ISO No. 514 |
| <span style="color: blue;">■</span> Blauer Ring/Blue ring/Bague bleue     | Standard/Standard/Standard   | 100 µm | ISO No. 524 |

# Diamantscheiben • Diamond Discs • Disques diamantés

Flex

806 104 ...

345



|                                |           |             |
|--------------------------------|-----------|-------------|
| Größe/ Size/ Taille            | ø 1/10 mm | 220         |
| Belegung/ Coating/ Granulation | B mm      | voll - full |
| Kopf/ Head Length/ Longueur mm | L mm      | 0,30        |

Shank   1


HP **806.104. ...** **345.524.220**

opt. 25.000 rpm

Flex TURBO

806 104 ...

362



|                                |           |             |             |
|--------------------------------|-----------|-------------|-------------|
| Größe/ Size/ Taille            | ø 1/10 mm | 080         | 100         |
| Belegung/ Coating/ Granulation | B mm      | voll - full | voll - full |
| Kopf/ Head Length/ Longueur mm | L mm      | 0,30        | 0,30        |

Shank   1

HP **TURBO 806.104. ...** **362.524.080** **362.524.100**

opt. 30.000 rpm



806.104.362.524.100

**Turbo Diamantscheibe** für die Laborturbine zum Separieren und Trennen von keramischen Werkstoffen.

**Turbo Diamond Disc** for use in a laboratory turbine for contouring interdental spaces and cutting through ceramics.

**Disque diamanté Turbo** pour la turbine de laboratoire servant à la finition des espaces inter-dentaires et à la séparation des matériaux céramiques.

Flex

806 314 ...

C8-FG



|                                |           |             |
|--------------------------------|-----------|-------------|
| Größe/ Size/ Taille            | ø 1/10 mm | 080         |
| Belegung/ Coating/ Granulation | B mm      | voll - full |
| Kopf/ Head Length/ Longueur mm | L mm      | 0,30        |

Shank   1

FG **TURBO 806.314.362.524.080** **C8-FG**

opt. 200.000 rpm



**FG**

C8-FG

|                                                                                     |                                     |                              |        |             |
|-------------------------------------------------------------------------------------|-------------------------------------|------------------------------|--------|-------------|
|  | Gelber Ring/Yellow Ring/Bague jaune | Superfein/Superfine/Superfin | 30 µm  | ISO No. 504 |
|  | Roter Ring/Red ring/Bague rouge     | Fein/Fine/Fin                | 50 µm  | ISO No. 514 |
|  | Blauer Ring/Blue ring/Bague bleue   | Standard/Standard/Standard   | 100 µm | ISO No. 524 |