





# Diamantscheiben • Diamond Discs • Disques diamantés

Superflex  
806 104 ...

351



|                                |           |             |
|--------------------------------|-----------|-------------|
| Größe/ Size/ Taille            | ø 1/10 mm | 220         |
| Belegung/ Coating/ Granulation | B mm      | voll - full |
| Kopf/ Head Length/ Longueur mm | L mm      | 0,10        |

|           |             |
|-----------|-------------|
| ø 1/10 mm | 220         |
| B mm      | voll - full |
| L mm      | 0,15        |

Shank   

HP 806.104. ...





351.514.220

351.524.220

opt. 25.000 rpm

Superflex  
806 104 ...

352



|                                |           |             |
|--------------------------------|-----------|-------------|
| Größe/ Size/ Taille            | ø 1/10 mm | 220         |
| Belegung/ Coating/ Granulation | B mm      | voll - full |
| Kopf/ Head Length/ Longueur mm | L mm      | 0,10        |

|           |             |
|-----------|-------------|
| ø 1/10 mm | 220         |
| B mm      | voll - full |
| L mm      | 0,15        |

Shank   

HP 806.104. ...





352.514.220

352.524.220

opt. 25.000 rpm

Superflex  
806 104 ...

353



|                                |             |             |             |
|--------------------------------|-------------|-------------|-------------|
| Größe/ Size/ Taille            | 220         | 220         | 220         |
| Belegung/ Coating/ Granulation | voll - full | voll - full | voll - full |
| Kopf/ Head Length/ Longueur mm | 0,12        | 0,15        | 0,25        |

Shank   

HP 806.104. ...

353.504.220

353.514.220

353.524.220

opt. 25.000 rpm



806.104.353.514.220

|                                                                                     |                                     |                              |        |             |
|-------------------------------------------------------------------------------------|-------------------------------------|------------------------------|--------|-------------|
|  | Gelber Ring/Yellow Ring/Bague jaune | Superfein/Superfine/Superfin | 30 µm  | ISO No. 504 |
|  | Roter Ring/Red ring/Bague rouge     | Fein/Fine/Fin                | 50 µm  | ISO No. 514 |
|  | Blauer Ring/Blue ring/Bague bleue   | Standard/Standard/Standard   | 100 µm | ISO No. 524 |